




NEW additions to our line of Refractive Drapes

MK-D360 Dual Refractive Drape

Our most versatile refractive drape



16" x 9.75", aperture area 7" x 3"




Ref. Number
MK-D360

Description
Dual Refractive Drape

Bag Qty.
20 per bag

MK-D1022 Refractive Drape

Split aperture for easy retraction of eye lid & lashes



16" x 10.7", aperture 1.5" x 1"

Ref. Number
MK-D1022

Description
Refractive Drape

Bag Qty.
20 per bag

From:


MICROspecialties, inc.
specialty ophthalmic products
www.microspecialties.com